

Data Path : S:\HPCHEM1\HPLC_N\Data\HN052313\
 Data File : HN005269.D
 Signal(s) : VWD1A.CH
 Acq On : 23-May-2013, 19:05
 Operator : QM
 Sample : E2136-01MS
 Misc :
 ALS Vial : 11 Sample Multiplier: 1

Instrument :
 HPLC_N
 ClientSampleId :
 GF004MS

Integration File: autoint1.e
 Quant Time: May 28 08:13:26 2013
 Quant Method : Z:\HPCHEM1\HPLC_N\METHOD\HN052213.M
 Quant Title :
 QLast Update : Fri May 24 05:23:33 2013
 Response via : Initial Calibration
 Integrator: ChemStation

Volume Inj. : 50 ul
 Signal Phase : ZORBAX Extend-C18
 Signal Info : 4.6 X 100mm 3.5-Micron 80A

	Compound	R.T.	Response	Conc Units

System Monitoring Compounds				
1) S	4-Nitroaniline	1.356	864514	203.403 ug/L
Target Compounds				
2) T	HMX	0.915	597206	140.376 ug/L
3) T	RDX	1.760	711520	138.334 ug/L
4) T	1,3,5-Trinitrobenzene	2.707	1569355	143.111 ug/L
5) T	1,3-Dinitrobenzene	3.596	2246288	172.486 ng/ml
6) T	3,5-Dinitroaniline	4.021	1840776	118.244 ug/L
7) T	Nitrobenzene	4.785	1334255	143.348 ug/L
8) T	Tetryl	6.455	621627	103.497 ug/L m
9) T	2,4,6-Trinitrotoluene	6.608	1806346	154.519 ug/L m
10) T	2-Amino-4,6-dinitroto...	7.259	1433704	134.276 ug/L
11) T	4-Amino-2,6-dinitroto...	7.623	948681	155.355 ug/L
12) T	2,4-Dinitrotoluene	8.273	2056950	141.410 ug/L
13) T	2,6-Dinitrotoluene	8.690	1186948	143.870 ug/L
14) T	2-Nitrotoluene	11.873	930836	149.248 ug/L m
15) T	4-Nitrotoluene	12.626	896260	170.865 ug/L
16) T	3-Nitrotoluene	13.947	1018470	153.344 ug/L
17) T	PETN	19.371	1294018	150.649 ug/L m

(f)=RT Delta > 1/2 Window

(m)=manual int.

Data Path : S:\HPCHEM1\HPLC_N\Data\HN052313\
Data File : HN005269.D
Signal(s) : VWD1A.CH
Acq On : 23-May-2013, 19:05
Operator : QM
Sample : E2136-01MS
Misc :
ALS Vial : 11 Sample Multiplier: 1

Instrument :
HPLC_N
ClientSampleId :
GF004MS

Integration File: autoint1.e
Quant Time: May 28 08:13:26 2013
Quant Method : Z:\HPCHEM1\HPLC_N\METHOD\HN052213.M
Quant Title :
QLast Update : Fri May 24 05:23:33 2013
Response via : Initial Calibration
Integrator: ChemStation

Volume Inj. : 50 ul
Signal Phase : ZORBAX Extend-C18
Signal Info : 4.6 X 100mm 3.5-Micron 80A

