

Method : H:\VOL1\HPLC\HPLC N\METHODS\HN101012.M (Chemstation Integrator)
 Title : SW8330 / 8330A confirmation
 Last Update : Thu Oct 11 16:07:34 2012

Calibration Files

CAL1 =HN000529.D CAL2 =HN000528.D CAL3 =HN000527.D
 CAL4 =HN000560.D CAL5 =HN000525.D CAL6 =HN000524.D

	Compound	CAL1	CAL2	CAL3	CAL4	CAL5	CAL6	Avg	%RSD
1) S	4-Nitroaniline	1.146	0.990	0.998	1.042	1.040	0.992	1.035 E4	5.75
2) T	HMX	1.050	1.048	1.018	1.033	1.027	1.019	1.033 E4	1.34
3) T	RDX	1.599	1.326	1.342	1.373	1.365	1.356	1.393 E4	7.34
4) T	1,3,5-Trinitrobenze	2.720	2.681	2.720	2.899	2.873	2.828	2.787 E4	3.29
5) T	Tetryl	2.356	2.100	2.057	2.153	2.149	2.151	2.161 E4	4.76
6) T	2,4,6-Trinitrotolue	2.787	2.538	2.772	2.685	2.679	2.650	2.685 E4	3.37
7) T	2-Amino-4,6-dinitro	4.371	4.672	4.681	4.810	4.773	4.704	4.669 E4	3.33
8) T	4-Amino-2,6-dinitro	4.371	4.672	4.681	4.810	4.773	4.704	4.669 E4	3.33
9) T	2,4-Dinitrotoluene	3.822	3.692	3.756	3.844	3.829	3.784	3.788 E4	1.51
10) T	2,6-Dinitrotoluene	1.881	2.029	1.934	1.991	1.978	1.955	1.961 E4	2.59
11) T	1,3-Dinitrobenzene	3.929	4.022	4.029	4.141	4.101	4.053	4.046 E4	1.80
12) T	Nitrobenzene	2.733	2.576	2.550	2.611	2.547	2.655	2.612 E4	2.75
13) T	2-Nitrotoluene	1.361	1.550	1.629	1.697	1.656	1.697	1.598 E4	8.03
14) T	3-Nitrotoluene	2.712	2.247	1.864	1.872	1.790	1.830	2.053 E4	17.69
15) T	4-Nitrotoluene	1.013	1.467	1.404	1.455	1.418	1.431	1.365 E4	12.75
16) T	3,5-Dinitroaniline	3.478	3.304	3.282	3.359	3.338	3.305	3.344 E4	2.12
17) T	Nitroqlvcerin		0.786	1.447	1.399	1.447	1.410	1.298 E2	22.08
18) T	PETN						0.000	0	-1.00