

Method : H:\VOL1\HPLC\HPLC N\METHODS\HN112012.M (Chemstation Integrator)
 Title : SW8330 / 8330A confirmation
 Last Update : Wed Nov 21 16:37:45 2012

Calibration Files

CAL1 =HN000595.D CAL2 =HN000594.D CAL3 =HN000593.D
 CAL4 =HN000592.D CAL5 =HN000591.D CAL6 =HN000589.D

	Compound	CAL1	CAL2	CAL3	CAL4	CAL5	CAL6	Avg	%RSD
1) S	4-Nitroaniline	1.090	0.981	0.994	1.003	0.983	0.998	1.008 E4	4.04
2) T	HMX	1.001	0.990	1.009	1.027	1.005	1.019	1.009 E4	1.30
3) T	RDX	1.479	1.340	1.375	1.364	1.349	1.367	1.379 E4	3.66
4) T	1,3,5-Trinitrobenze	2.743	2.780	2.807	2.826	2.788	2.815	2.793 E4	1.07
5) T	Tetryl	2.170	2.210	2.148	2.157	2.128	2.129	2.157 E4	1.42
6) T	2,4,6-Trinitrotolue	2.822	2.828	2.668	2.642	2.642	2.655	2.710 E4	3.32
7) T	2-Amino-4,6-dinitro	4.752	4.916	4.793	4.763	4.691	4.700	4.769 E4	1.71
8) T	4-Amino-2,6-dinitro	4.919	4.916	4.793	4.763	4.691	4.700	4.797 E4	2.10
9) T	2,4-Dinitrotoluene	4.048	3.871	3.793	3.816	3.739	3.783	3.842 E4	2.86
10) T	2,6-Dinitrotoluene	2.306	2.067	1.954	1.964	1.936	1.964	2.032 E4	6.98
11) T	1,3-Dinitrobenzene	4.730	4.370	4.096	4.076	4.030	4.047	4.225 E4	6.57
12) T	Nitrobenzene	2.683	2.547	2.600	2.653	2.602	2.641	2.621 E4	1.84
13) T	2-Nitrotoluene	1.562	1.772	1.724	1.708	1.686	1.690	1.690 E4	4.15
14) T	3-Nitrotoluene	1.719	1.918	1.827	1.825	1.784	1.801	1.812 E4	3.58
15) T	4-Nitrotoluene	1.142	1.493	1.455	1.439	1.417	1.419	1.394 E4	9.09
16) T	3,5-Dinitroaniline	3.579	3.265	3.294	3.332	3.257	3.300	3.338 E4	3.63
17) T	Nitroqlvcerin		1.478	1.354	1.439	1.379	1.363	1.403 E2	3.81
18) T	PETN						0.000	0	-1.00