

Data File : Z:\HPCHEM1\HPLC\_N\DATA\HN052914\HN006053.D Vial: 2  
Acq On : 29-May-2014, 14:09 Operator: HNP  
Sample : HSTDICC20 Inst : G1314A  
Misc : Multiplr: 1.00  
IntFile : AUTOINT1.E  
Quant Time: May 30 3:35 2014 Quant Results File: HN0529NG.RES

Instrument :  
HPLC\_N  
LabSampleId :  
HSTDICC20

Quant Method : Z:\HPCHEM1\HPLC\_N\METHODS\HN0529NG.M  
Title :  
Last Update : Fri May 30 03:32:30 2014  
Response via : Initial Calibration  
DataAcq Meth : 8330ANG.M

Volume Inj. : 50 ul  
Signal Phase : ZORBAX Extend-C18  
Signal Info : 4.6 X 100mm 3.5-Micron 80A

| Compound                    | R.T. | Response  | Conc Units   |
|-----------------------------|------|-----------|--------------|
| -----                       |      |           |              |
| System Monitoring Compounds |      |           |              |
| 1) S 4-Nitroaniline         | 1.42 | 74237351  | 21.059 ng/ul |
| Target Compounds            |      |           |              |
| 2) T Nitroglycerin          | 6.58 | 144999711 | 19.145 ng/ul |

-----  
Analyst Signature:\_\_\_\_\_ Analyst Name:\_\_\_\_\_ Date:\_\_\_\_\_

-----REASONS FOR MANUAL INTEGRATIONS-----

\_\_\_Poor resolution of peaks exhibited on chromatogram.Compound #:\_\_\_\_\_

\_\_\_Peak Integrated by software incorrectly.Compound #:\_\_\_\_\_

\_\_\_OTHER:\_\_\_\_\_Compound #:\_\_\_\_\_

(f)=RT Delta > 1/2 Window (m)=manual int.

HN006053.D HN0529NG.M

Wed Jul 16 17:54:24 2014

RPT1

# Quantitation Report

Data File : Z:\HPCHEM1\HPLC\_N\DATA\HN052914\HN006053.D Vial: 2  
Acq On : 29-May-2014, 14:09 Operator: HNP  
Sample : HSTDICC20 Inst : G1314A  
Misc : Multiplr: 1.00  
IntFile : AUTOINT1.E  
Quant Time: May 30 3:35 2014 Quant Results File: HN0529NG.RES

Instrument :  
HPLC\_N  
LabSampleId :  
HSTDICC20

Quant Method : Z:\HPCHEM1\HPLC\_N\METHODS\HN0529NG.M  
Title :  
Last Update : Fri May 30 03:32:30 2014  
Response via : Multiple Level Calibration  
DataAcq Meth : 8330ANG.M

Volume Inj. : 50 ul  
Signal Phase : ZORBAX Extend-C18  
Signal Info : 4.6 X 100mm 3.5-Micron 80A

