

Method : H:\VOL1\HPLC\HPLC N\METHODS\HN081512.M (Chemstation Integrator)
 Title : SW8330 / 8330A confirmation
 Last Update : Thu Aug 16 10:32:31 2012

Calibration Files

CAL1 =HN000452.D CAL2 =HN000451.D CAL3 =HN000450.D
 CAL4 =HN000449.D CAL5 =HN000448.D CAL6 =HN000447.D

	Compound	CAL1	CAL2	CAL3	CAL4	CAL5	CAL6	Avg	%RSD
1) S	4-Nitroaniline	0.954	0.982	0.975	1.018	0.991	0.985	0.984 E4	2.12
2) T	HMX	1.064	0.989	0.978	1.068	1.013	1.026	1.023 E4	3.66
3) T	RDX	1.311	1.332	1.329	1.405	1.351	1.351	1.346 E4	2.41
4) T	1,3,5-Trinitrobenze	2.589	3.447	2.890	2.910	2.850	2.828	2.919 E4	9.71
5) T	Tetryl	2.153	2.160	2.110	2.235	2.176	2.173	2.168 E4	1.88
6) T	2,4,6-Trinitrotolue	2.640	2.536	2.701	2.794	2.705	2.670	2.674 E4	3.19
7) T	2-Amino-4,6-dinitro	4.410	4.744	4.726	5.003	4.828	4.822	4.756 E4	4.11
8) T	4-Amino-2,6-dinitro	4.410	4.744	4.726	5.003	4.828	4.822	4.756 E4	4.11
9) T	2,4-Dinitrotoluene	3.686	3.739	3.720	3.923	3.779	3.771	3.770 E4	2.19
10) T	2,6-Dinitrotoluene	1.883	1.983	1.922	2.027	1.952	1.947	1.953 E4	2.54
11) T	1,3-Dinitrobenzene	3.582	3.736	3.787	4.015	3.898	3.902	3.820 E4	3.98
12) T	Nitrobenzene	2.649	2.529	2.547	2.714	2.620	2.634	2.616 E4	2.62
13) T	2-Nitrotoluene	1.510	1.638	1.640	1.738	1.681	1.683	1.648 E4	4.67
14) T	3-Nitrotoluene	1.781	1.877	1.795	1.878	1.812	1.815	1.826 E4	2.27
15) T	4-Nitrotoluene	1.590	1.470	1.414	1.488	1.435	1.435	1.472 E4	4.33
16) T	3,5-Dinitroaniline	3.491	3.302	3.246	3.451	3.295	3.326	3.352 E4	2.88
17) T	Nitroqllycerin							0.000 0	-1.00
18) T	PETN							0.000 0	-1.00