

Method : H:\VOL1\HPLC\HPLC N\METHODS\HN081512.M (Chemstation Integrator)
 Title : SW8330 / 8330A confirmation
 Last Update : Fri Aug 17 11:49:03 2012
 Response via : Initial Calibration

Continuing Calibration File: HN000453.D

Min. RRF : 0.000 Min. Rel. Area : 50%
 Max. RRF Dev : 20% Max. Rel. Area : 150%

	Compound	AvgRF	CCRF		%Dev	Area%
1 S	4-Nitroaniline	9.844	0.000	E3	100.0#	-0#
2 T	HMX	10.230	10.341	E3	-1.1	97
3 T	RDX	13.464	14.149	E3	-5.1	101
4 T	1,3,5-Trinitrobenzene	29.189	29.829	E3	-2.2	205#
5 T	Tetryl	21.677	22.301	E3	-2.9	100
6 T	2,4,6-Trinitrotoluene	26.744	28.488	E3	-6.5	204#
7 T	2-Amino-4,6-dinitrotoluene	47.555	49.603	E3	-4.3	198#
8 T	4-Amino-2,6-dinitrotoluene	47.555	49.603	E3	-4.3	198#
9 T	2,4-Dinitrotoluene	37.699	37.235	E3	1.2	190#
10 T	2,6-Dinitrotoluene	19.525	19.815	E3	-1.5	195#
11 T	1,3-Dinitrobenzene	38.201	39.497	E3	-3.4	197#
12 T	Nitrobenzene	26.156	26.523	E3	-1.4	195#
13 T	2-Nitrotoluene	16.484	17.228	E3	-4.5	99
14 T	3-Nitrotoluene	18.265	18.327	E3	-0.3	98
15 T	4-Nitrotoluene	14.720	14.642	E3	0.5	98
16 T	3,5-Dinitroaniline	33.519	33.080	E3	1.3	96
17 T	Nitroglycerin	0.000	0.000		0.0	0#
18 T	PETN	0.000	0.000		0.0	0#