

Signal #1 : S:\HPCHEM1\HPLC_N\DATA\HN093014\HN006241.D\VWD1A.CH Vial: Instrument :
Signal #2 : S:\HPCHEM1\HPLC_N\DATA\HN093014\HN006241.D\VWD1A.CH HPLC_N
Acq On : 30-Sep-2014, 11:50 Operator: HNP ClientSampleId :
Sample : HSTDICC2000 Inst : G1314A HSTDICC2000
Misc : Multiplr: 1.00
IntFile Signal #1: autoint1.e IntFile Signal #2: autoint2.e
Quant Time: Oct 1 5:55 2014 Quant Results File: HN093014.RES

Quant Method : Z:\HPCHEM1\HPLC_N\METHODS\HN093014.M
Title :
Last Update : Wed Oct 01 04:09:40 2014
Response via : Initial Calibration
DataAcq Meth : 8330AGIP.M

Volume Inj. : 50 ul
Signal #1 Phase : ZORBAX Extend-C18 Signal #2 Phase:
Signal #1 Info : 4.6 X 100mm 3.5-M Signal #2 Info :

Compound	R.T.	Response	Conc Units

System Monitoring Compounds			
1) S 4-Nitroaniline	1.07	7010204	1149.950 ug/L
Target Compounds			
2) T HMX	0.50	2572812	291.657 ug/L
3) T RDX	2.15	8671225	723.139 ug/L
4) T 1,3,5-Trinitrobenzene	3.54	17321127	1231.274 ug/L
5) T 1,3-Dinitrobenzene	4.75	26666710	1237.214 ng/ml
6) T 3,5-Dinitroaniline	5.27	20038902	1185.708 ug/L
7) T Nitrobenzene	6.46	16853140	1275.379 ug/L

Analyst Signature:_____ Analyst Name:_____ Date:_____

-----REASONS FOR MANUAL INTEGRATIONS-----

___Poor resolution of peaks exhibited on chromatogram.Compound #:_____

___Peak Integrated by software incorrectly.Compound #:_____

___OTHER:_____Compound #:_____

(f)=RT Delta > 1/2 Window (m)=manual int.

HN006241.D HN093014.M

Thu Oct 09 13:56:15 2014

RPT1

Quantitation Report

Signal #1 : S:\HPCHEM1\HPLC_N\DATA\HN093014\HN006241.D\VWD1A.CH Vial:
Signal #2 : S:\HPCHEM1\HPLC_N\DATA\HN093014\HN006241.D\VWD1A.CH
Acq On : 30-Sep-2014, 11:50 Operator: HNP
Sample : HSTDICC2000 Inst : G1314A
Misc : Multiplr: 1.00
IntFile Signal #1: autoint1.e IntFile Signal #2: autoint2.e
Quant Time: Oct 1 5:55 2014 Quant Results File: HN093014.RES

Instrument :
HPLC_N
ClientSampleId :
HSTDICC2000

Quant Method : Z:\HPCHEM1\HPLC_N\METHODS\HN093014.M
Title :
Last Update : Wed Oct 01 04:09:40 2014
Response via : Multiple Level Calibration
DataAcq Meth : 8330AGIP.M

Volume Inj. : 50 ul
Signal #1 Phase : ZORBAX Extend-C18 Signal #2 Phase:
Signal #1 Info : 4.6 X 100mm 3.5-M Signal #2 Info :

