

Signal #1 : Z:\HPCHEM1\HPLC_N\DATA\HN101514\HN006313.D\VWD1A.CH Vial:
 Signal #2 : Z:\HPCHEM1\HPLC_N\DATA\HN101514\HN006313.D\VWD1A.CH
 Acq On : 15-Oct-2014, 17:31 Operator: HNP
 Sample : F4222-12 Inst : G1314A
 Misc : Multiplr: 1.00
 IntFile Signal #1: autoint1.e IntFile Signal #2: autoint2.e
 Quant Time: Oct 16 11:05 2014 Quant Results File: HN100214.RES

Instrument :
 HPLC_N
 ClientSampleId :
 YS19-EB02-100814-GW

Manual Integrations
 APPROVED

mohammad
 10/21/2014 10:20:11 AM

Quant Method : Z:\HPCHEM1\HPLC_N\METHODS\HN100214.M
 Title :
 Last Update : Fri Oct 03 10:34:01 2014
 Response via : Initial Calibration
 DataAcq Meth : 8330AGIP.M

Volume Inj. : 50 ul
 Signal #1 Phase : ZORBAX Extend-C18 Signal #2 Phase:
 Signal #1 Info : 4.6 X 100mm 3.5-M Signal #2 Info :

Compound	R.T.	Response	Conc Units

System Monitoring Compounds			
1) S 4-Nitroaniline	1.31	1146321	288.417 ug/L m

Target Compounds

 Analyst Signature:_____ Analyst Name:_____ Date:_____

-----REASONS FOR MANUAL INTEGRATIONS-----

___Poor resolution of peaks exhibited on chromatogram.Compound #:_____

___Peak Integrated by software incorrectly.Compound #:_____

___OTHER:_____Compound #:_____

(f)=RT Delta > 1/2 Window (m)=manual int.

HN006313.D HN100214.M Fri Oct 17 18:56:40 2014 RPT1

Quantitation Report

Signal #1 : Z:\HPCHEM1\HPLC_N\DATA\HN101514\HN006313.D\VWD1A.CH Vial:
Signal #2 : Z:\HPCHEM1\HPLC_N\DATA\HN101514\HN006313.D\VWD1A.CH
Acq On : 15-Oct-2014, 17:31 Operator: HNP
Sample : F4222-12 Inst : G1314A
Misc : Multiplr: 1.00
IntFile Signal #1: autoint1.e IntFile Signal #2: autoint2.e
Quant Time: Oct 16 11:05 2014 Quant Results File: HN100214.RES

Instrument :
HPLC_N
ClientSampleId :
YS19-EB02-100814-GW

Manual Integrations
APPROVED

mohammad
10/21/2014 10:20:11 AM

Quant Method : Z:\HPCHEM1\HPLC_N\METHODS\HN100214.M
Title :
Last Update : Fri Oct 03 10:34:01 2014
Response via : Multiple Level Calibration
DataAcq Meth : 8330AGIP.M

Volume Inj. : 50 ul
Signal #1 Phase : ZORBAX Extend-C18 Signal #2 Phase:
Signal #1 Info : 4.6 X 100mm 3.5-M Signal #2 Info :

