

Signal #1 : X:\HPCHEM1\HPLC N\DATA\HN101514\HN006327.D\VWD1A.CH Vial: 26
 Signal #2 : X:\HPCHEM1\HPLC N\DATA\HN101514\HN006327.D\VWD1A.CH
 Acq On : 16-Oct-2014, 00:03 Operator: HNP
 Sample : F4244-10 Inst : G1314A
 Misc : Multiplr: 1.00
 IntFile Signal #1: autoint1.e IntFile Signal #2: autoint2.e
 Quant Time: Oct 16 11:41 2014 Quant Results File: HN100214.RES

Instrument :
 HPLC_N
 ClientSampleId :
 YS19-GW12-1014A

Manual Integrations
 APPROVED

mohammad
 10/21/2014 10:20:42 AM

Quant Method : Z:\HPCHEM1\HPLC_N\METHODS\HN100214.M
 Title :
 Last Update : Fri Oct 03 10:34:01 2014
 Response via : Initial Calibration
 DataAcq Meth : 8330AGIP.M

Volume Inj. : 50 ul
 Signal #1 Phase : ZORBAX Extend-C18 Signal #2 Phase:
 Signal #1 Info : 4.6 X 100mm 3.5-M Signal #2 Info :

Compound	R.T.	Response	Conc Units

System Monitoring Compounds			
1) S 4-Nitroaniline	1.31	1040905	261.894 ug/L m
Target Compounds			

 Analyst Signature: _____ Analyst Name: _____ Date: _____

-----REASONS FOR MANUAL INTEGRATIONS-----

Poor resolution of peaks exhibited on chromatogram.Compound #:

Peak Integrated by software incorrectly.Compound #:

OTHER: _____ Compound #: _____

(f)=RT Delta > 1/2 Window (m)=manual int.

HN006327.D HN100214.M Thu Oct 30 05:15:11 2014

Quantitation Report

Signal #1 : X:\HPCHEM1\HPLC N\DATA\HN101514\HN006327.D\VWD1A.CH Vial: 26
 Signal #2 : X:\HPCHEM1\HPLC N\DATA\HN101514\HN006327.D\VWD1A.CH
 Acq On : 16-Oct-2014, 00:03 Operator: HNP
 Sample : F4244-10 Inst : G1314A
 Misc : Multiplr: 1.00
 IntFile Signal #1: autoint1.e IntFile Signal #2: autoint2.e
 Quant Time: Oct 16 11:41 2014 Quant Results File: HN100214.RES

Instrument :
 HPLC_N
 ClientSampleId :
 YS19-GW12-1014A

Manual Integrations
 APPROVED

mohammad
 10/21/2014 10:20:42 AM

Quant Method : Z:\HPCHEM1\HPLC_N\METHODS\HN100214.M
 Title :
 Last Update : Fri Oct 03 10:34:01 2014
 Response via : Multiple Level Calibration
 DataAcq Meth : 8330AGIP.M

Volume Inj. : 50 ul
 Signal #1 Phase : ZORBAX Extend-C18 Signal #2 Phase:
 Signal #1 Info : 4.6 X 100mm 3.5-M Signal #2 Info :

