

Signal #1 : X:\HPCHEM1\HPLC N\DATA\HN101714\HN006469.D\VWD1A.CH Vial: 44
 Signal #2 : X:\HPCHEM1\HPLC N\DATA\HN101714\HN006469.D\VWD1A.CH
 Acq On : 18-Oct-2014, 19:13 Operator: HNP
 Sample : F4322-17 Inst : G1314A
 Misc : Multiplr: 1.00
 IntFile Signal #1: AUTOINT1.E IntFile Signal #2: autoint2.e
 Quant Time: Oct 22 12:45 2014 Quant Results File: HN1016NG.RES

Instrument :
 HPLC_N
 ClientSampleId :
 YS19-EB01-101514-SO

Manual Integrations
 APPROVED

mohammad
 10/30/2014 10:27:33 AM

Quant Method : Z:\HPCHEM1\HPLC_N\METHODS\HN1016NG.M
 Title :
 Last Update : Fri Oct 17 04:15:53 2014
 Response via : Initial Calibration
 DataAcq Meth : 8330ANG.M

Volume Inj. : 50 ul
 Signal #1 Phase : ZORBAX Extend-C18 Signal #2 Phase:
 Signal #1 Info : 4.6 X 100mm 3.5-M Signal #2 Info :

Compound	R.T.	Response	Conc Units

System Monitoring Compounds			
1) S 4-Nitroaniline	1.31	1216568	0.306 ng/ulm

Target Compounds

 Analyst Signature: _____ Analyst Name: _____ Date: _____

-----REASONS FOR MANUAL INTEGRATIONS-----

Poor resolution of peaks exhibited on chromatogram.Compound #:

Peak Integrated by software incorrectly.Compound #:

OTHER:

Compound #: _____

(f)=RT Delta > 1/2 Window

(m)=manual int.

HN006469.D HN1016NG.M

Thu Oct 30 10:05:34 2014

Quantitation Report

Signal #1 : X:\HPCHEM1\HPLC N\DATA\HN101714\HN006469.D\VWD1A.CH Vial: 44
 Signal #2 : X:\HPCHEM1\HPLC N\DATA\HN101714\HN006469.D\VWD1A.CH
 Acq On : 18-Oct-2014, 19:13 Operator: HNP
 Sample : F4322-17 Inst : G1314A
 Misc : Multiplr: 1.00
 IntFile Signal #1: AUTOINT1.E IntFile Signal #2: autoint2.e
 Quant Time: Oct 22 12:45 2014 Quant Results File: HN1016NG.RES

Instrument :
 HPLC_N
 ClientSampleId :
 YS19-EB01-101514-SO

Manual Integrations
 APPROVED

mohammad
 10/30/2014 10:27:33 AM

Quant Method : Z:\HPCHEM1\HPLC_N\METHODS\HN1016NG.M
 Title :
 Last Update : Fri Oct 17 04:15:53 2014
 Response via : Multiple Level Calibration
 DataAcq Meth : 8330ANG.M

Volume Inj. : 50 ul
 Signal #1 Phase : ZORBAX Extend-C18 Signal #2 Phase:
 Signal #1 Info : 4.6 X 100mm 3.5-M Signal #2 Info :

