

Signal #1 : X:\HPCHEM1\HPLC N\DATA\HN101714\HN006470.D\VWD1A.CH Vial: 45
 Signal #2 : X:\HPCHEM1\HPLC N\DATA\HN101714\HN006470.D\VWD1A.CH
 Acq On : 18-Oct-2014, 19:25 Operator: HNP
 Sample : F4322-18 Inst : G1314A
 Misc : Multiplr: 1.00
 IntFile Signal #1: AUTOINT1.E IntFile Signal #2: autoint2.e
 Quant Time: Oct 22 12:47 2014 Quant Results File: HN1016NG.RES

Instrument :
 HPLC_N
 ClientSampleId :
 YS19-EB02-101514-SO

Manual Integrations
 APPROVED

mohammad
 10/30/2014 10:27:35 AM

Quant Method : Z:\HPCHEM1\HPLC_N\METHODS\HN1016NG.M
 Title :
 Last Update : Fri Oct 17 04:15:53 2014
 Response via : Initial Calibration
 DataAcq Meth : 8330ANG.M

Volume Inj. : 50 ul
 Signal #1 Phase : ZORBAX Extend-C18 Signal #2 Phase:
 Signal #1 Info : 4.6 X 100mm 3.5-M Signal #2 Info :

Compound	R.T.	Response	Conc Units

System Monitoring Compounds			
1) S 4-Nitroaniline	1.31	1272023	0.320 ng/ulm

Target Compounds

 Analyst Signature: _____ Analyst Name: _____ Date: _____

-----REASONS FOR MANUAL INTEGRATIONS-----

Poor resolution of peaks exhibited on chromatogram.Compound #:

Peak Integrated by software incorrectly.Compound #:

OTHER:

Compound #: _____

(f)=RT Delta > 1/2 Window

(m)=manual int.

HN006470.D HN1016NG.M

Thu Oct 30 10:05:40 2014

Quantitation Report

Signal #1 : X:\HPCHEM1\HPLC N\DATA\HN101714\HN006470.D\VWD1A.CH Vial: 45

Signal #2 : X:\HPCHEM1\HPLC N\DATA\HN101714\HN006470.D\VWD1A.CH

Acq On : 18-Oct-2014, 19:25

Operator: HNP

Sample : F4322-18

Inst : G1314A

Misc :

Multiplr: 1.00

IntFile Signal #1: AUTOINT1.E

IntFile Signal #2: autoint2.e

Quant Time: Oct 22 12:47 2014 Quant Results File: HN1016NG.RES

Instrument :

HPLC_N

ClientSampleId :

YS19-EB02-101514-SO

Manual Integrations
APPROVEDmohammad
10/30/2014 10:27:35 AM

Quant Method : Z:\HPCHEM1\HPLC_N\METHODS\HN1016NG.M

Title :

Last Update : Fri Oct 17 04:15:53 2014

Response via : Multiple Level Calibration

DataAcq Meth : 8330ANG.M

Volume Inj. : 50 ul

Signal #1 Phase : ZORBAX Extend-C18 Signal #2 Phase:

Signal #1 Info : 4.6 X 100mm 3.5-M Signal #2 Info :

